

Y-box-binding protein 3

Catalog No: #AP88944



Package Size: #AP88944-1 50ug #AP88944-2 100ug #AP88944-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Y-box-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cold shock domain-containing protein A,DNA-binding protein A,Y-box protein 3,Csda,Msy4
Accession No.	Uniprot ID:Q9JKB3
Uniprot	Q9JKB3
GeneID	56449;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ybx3

Note: This product is for in vitro research use only