

Ubiquitin-associated protein 2

Catalog No: #AP88945



Package Size: #AP88945-1 50ug #AP88945-2 100ug #AP88945-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	UBAP-2,Protein lingerer homolog 1,mLig-1,Lig1
Accession No.	Uniprot ID:Q91VX2
Uniprot	Q91VX2
GeneID	68926;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ubap2

Note: This product is for in vitro research use only