

## Ubiquitin-conjugating enzyme E2 A

Catalog No: #AP88947



Package Size: #AP88947-1 50ug #AP88947-2 100ug #AP88947-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubiquitin-conjugating enzyme E2 A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	E2 ubiquitin-conjugating enzyme A,RAD6 homolog A,HR6A,Ubiquitin carrier protein A,Ubiquitin-protein ligase A,Rad6a
Accession No.	Uniprot ID:Q9Z255
Uniprot	Q9Z255
GeneID	22209;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ube2a

Note: This product is for in vitro research use only