

# Lactotransferrin

Catalog No: #AP88949



Package Size: #AP88949-1 50ug #AP88949-2 100ug #AP88949-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lactotransferrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lactoferrin,Growth-inhibiting protein 12,Talalactoferrin,GIG12,LF
Accession No.	Uniprot ID:P02788
Uniprot	P02788
GeneID	4057;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LTF

Note: This product is for in vitro research use only