

Tripartite motif-containing protein 3

Catalog No: #AP88950



Package Size: #AP88950-1 50ug #AP88950-2 100ug #AP88950-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tripartite motif-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Brain-expressed RING finger protein,RING finger protein 22,Berp,Rnf22
Accession No.	Uniprot ID:O70277
Uniprot	O70277
GeneID	83616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trim3

Note: This product is for in vitro research use only