

Vascular endothelial growth factor A

Catalog No: #AP88963



Package Size: #AP88963-1 50ug #AP88963-2 100ug #AP88963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	VEGF-A,Vascular permeability factor,VPF,VEGF
Accession No.	Uniprot ID:Q9MYV3
Uniprot	Q9MYV3
GeneID	403802;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VEGFA

Note: This product is for in vitro research use only