

Zinc finger protein 335

Catalog No: #AP88964



Package Size: #AP88964-1 50ug #AP88964-2 100ug #AP88964-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein 335
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NRC-interacting factor 1,NIF-1,Zfp335
Accession No.	Uniprot ID:A2A5K6
Uniprot	A2A5K6
GeneID	329559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Znf335

Note: This product is for in vitro research use only