

## Tenascin-N

Catalog No: #AP88966



Package Size: #AP88966-1 50ug #AP88966-2 100ug #AP88966-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tenascin-N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TN-N,Tenascin-W,TN-W,TNW
Accession No.	Uniprot ID:Q9UQP3
Uniprot	Q9UQP3
GeneID	63923;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNN

Note: This product is for in vitro research use only