

Visinin-like protein 1

Catalog No: #AP88969



Package Size: #AP88969-1 50ug #AP88969-2 100ug #AP88969-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Visinin-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VILIP,Hippocalcin-like protein 3,HLP3,VISL1
Accession No.	Uniprot ID:P62760
Uniprot	P62760
GeneID	7447;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VSNL1

Note: This product is for in vitro research use only