

## L-dopachrome tautomerase

Catalog No: #AP88970



Package Size: #AP88970-1 50ug #AP88970-2 100ug #AP88970-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-dopachrome tautomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	DCT,L-dopachrome Delta-isomerase,Tyrosinase-related protein 2,TRP-2,TYRP2
Accession No.	Uniprot ID:Q4R1H1
Uniprot	Q4R1H1
GeneID	574066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:DCT

Note: This product is for in vitro research use only