

Thymidylate synthase

Catalog No: #AP88971



Package Size: #AP88971-1 50ug #AP88971-2 100ug #AP88971-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Thymidylate synthase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | TS |
| Accession No. | Uniprot ID:P07607 |
| Uniprot | P07607 |
| GeneID | 22171; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tyms

Note: This product is for in vitro research use only