

Upstream stimulatory factor 1

Catalog No: #AP88986



Package Size: #AP88986-1 50ug #AP88986-2 100ug #AP88986-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Upstream stimulatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class B basic helix-loop-helix protein 11,bHLHb11,Major late transcription factor 1,BHLHB11,USF
Accession No.	Uniprot ID:P22415
Uniprot	P22415
GeneID	7391;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:USF1

Note: This product is for in vitro research use only