

Neurosecretory protein VGF

Catalog No: #AP88988



Package Size: #AP88988-1 50ug #AP88988-2 100ug #AP88988-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurosecretory protein VGF
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VGf8a protein
Accession No.	Uniprot ID:P20156
Uniprot	P20156
GeneID	29461;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Vgf

Note: This product is for in vitro research use only