

T-box transcription factor TBX3

Catalog No: #AP88989



Package Size: #AP88989-1 50ug #AP88989-2 100ug #AP88989-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | T-box transcription factor TBX3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | T-box protein 3 |
| Accession No. | Uniprot ID:O73718 |
| Uniprot | O73718 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TBX3

Note: This product is for in vitro research use only