

Telomerase reverse transcriptase

Catalog No: #AP88990



Package Size: #AP88990-1 50ug #AP88990-2 100ug #AP88990-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Telomerase reverse transcriptase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Telomerase catalytic subunit
Accession No.	Uniprot ID:Q27ID4
Uniprot	Q27ID4
GeneID	518884;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TERT

Note: This product is for in vitro research use only