

Y-box-binding protein 3

Catalog No: #AP89002



Package Size: #AP89002-1 50ug #AP89002-2 100ug #AP89002-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Y-box-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cold shock domain-containing protein A,DNA-binding protein A,Single-strand DNA-binding protein NF-GMB,CSDA,DBPA
Accession No.	Uniprot ID:P16989
Uniprot	P16989
GeneID	8531;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:YBX3

Note: This product is for in vitro research use only