

GATA-type zinc finger protein 1

Catalog No: #AP89003



Package Size: #AP89003-1 50ug #AP89003-2 100ug #AP89003-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GATA-type zinc finger protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GATA-like protein 1, GLP-1, GLP1
Accession No.	Uniprot ID:P0C6A0
Uniprot	P0C6A0
GeneID	100125288;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ZGLP1

Note: This product is for in vitro research use only