

Vasoactive intestinal polypeptide receptor 2

Catalog No: #AP89004



Package Size: #AP89004-1 50ug #AP89004-2 100ug #AP89004-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vasoactive intestinal polypeptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VIP-R-2,Pituitary adenylate cyclase-activating polypeptide type III receptor,PACAP type III receptor
Accession No.	Uniprot ID:P35000
Uniprot	P35000
GeneID	29555;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vipr2

Note: This product is for in vitro research use only