

Trem-like transcript 2 protein

Catalog No: #AP89005



Package Size: #AP89005-1 50ug #AP89005-2 100ug #AP89005-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trem-like transcript 2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TLT-2, Triggering receptor expressed on myeloid cells-like protein 2, Tlt2
Accession No.	Uniprot ID:Q2LA85
Uniprot	Q2LA85
GeneID	328833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Trem12

Note: This product is for in vitro research use only