

Transferrin receptor protein 1

Catalog No: #AP89006



Package Size: #AP89006-1 50ug #AP89006-2 100ug #AP89006-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transferrin receptor protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TR,CD71,Trfr
Accession No.	Uniprot ID:Q99376
Uniprot	Q99376
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tfrc

Note: This product is for in vitro research use only