

# DNA topoisomerase 2-beta

Catalog No: #AP89010



Package Size: #AP89010-1 50ug #AP89010-2 100ug #AP89010-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA topoisomerase 2-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNA topoisomerase II, beta isozyme
Accession No.	Uniprot ID:Q02880
Uniprot	Q02880
GeneID	7155;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TOP2B

Note: This product is for in vitro research use only