

Tetraspanin-14

Catalog No: #AP89012



Package Size: #AP89012-1 50ug #AP89012-2 100ug #AP89012-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tetraspanin-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Tspan-14,DC-TM4F2,Transmembrane 4 superfamily member 14,TM4SF14
Accession No.	Uniprot ID:Q8NG11
Uniprot	Q8NG11
GeneID	81619;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TSPAN14

Note: This product is for in vitro research use only