

# Tumor necrosis factor receptor superfamily member 11A

Catalog No: #AP89016

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89016-1 50ug #AP89016-2 100ug #AP89016-3 1mg

## Description

Product Name	Tumor necrosis factor receptor superfamily member 11A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Osteoclast differentiation factor receptor,ODFR,Receptor activator of NF-KB,CD265,Rank
Accession No.	Uniprot ID:O35305
Uniprot	O35305
GeneID	21934;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tnfrsf11a

Note: This product is for in vitro research use only