

## Tumor necrosis factor receptor superfamily member 5

Catalog No: #AP89017



Package Size: #AP89017-1 50ug #AP89017-2 100ug #AP89017-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B-cell surface antigen CD40,Bp50,CD40L receptor,CDw40,CD40,TNFRSF5
Accession No.	Uniprot ID:P25942
Uniprot	P25942
GeneID	958;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD40

Note: This product is for in vitro research use only