

Triggering receptor expressed on myeloid cells 1

Catalog No: #AP89019



Package Size: #AP89019-1 50ug #AP89019-2 100ug #AP89019-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Triggering receptor expressed on myeloid cells 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TREM-1,CD354
Accession No.	Uniprot ID:Q6QUN5
Uniprot	Q6QUN5
GeneID	404547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TREM1

Note: This product is for in vitro research use only