

AN1-type zinc finger protein 6

Catalog No: #AP89020



Package Size: #AP89020-1 50ug #AP89020-2 100ug #AP89020-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AN1-type zinc finger protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Zinc finger A20 domain-containing protein 3,Za20d3
Accession No.	Uniprot ID:Q6DGF4
Uniprot	Q6DGF4
GeneID	293067;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zfand6

Note: This product is for in vitro research use only