

Zinc finger homeobox protein 4

Catalog No: #AP89021



Package Size: #AP89021-1 50ug #AP89021-2 100ug #AP89021-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Zinc finger homeobox protein 4 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Zinc finger homeodomain protein 4,ZFH-4,Zfh4 |
| Accession No. | Uniprot ID:Q9JJN2 |
| Uniprot | Q9JJN2 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zfhx4

Note: This product is for in vitro research use only