

Zinc finger protein 335

Catalog No: #AP89024



Package Size: #AP89024-1 50ug #AP89024-2 100ug #AP89024-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein 335
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NRC-interacting factor 1,NIF-1,Zfp335
Accession No.	Uniprot ID:G3V893
Uniprot	G3V893
GeneID	259270;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Znf335

Note: This product is for in vitro research use only