

# Sex hormone-binding globulin

Catalog No: #AP89025



Package Size: #AP89025-1 50ug #AP89025-2 100ug #AP89025-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sex hormone-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SHBG, Sex steroid-binding protein, SBP, Testis-specific androgen-binding protein, ABP
Accession No.	Uniprot ID:P08689
Uniprot	P08689
GeneID	24775;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Shbg

Note: This product is for in vitro research use only