

## Src-like-adapter 2

Catalog No: #AP89026



Package Size: #AP89026-1 50ug #AP89026-2 100ug #AP89026-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Src-like-adapter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Modulator of antigen receptor signaling,MARS,Src-like adapter protein 2,SLAP-2,C20orf156,SLAP2
Accession No.	Uniprot ID:Q9H6Q3
Uniprot	Q9H6Q3
GeneID	84174;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLA2

Note: This product is for in vitro research use only