

Vascular endothelial growth factor D

Catalog No: #AP89031



Package Size: #AP89031-1 50ug #AP89031-2 100ug #AP89031-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VEGF-D,c-Fos-induced growth factor,FIGF,Figf
Accession No.	Uniprot ID:O35251
Uniprot	O35251
GeneID	360457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vegfd

Note: This product is for in vitro research use only