

Transforming growth factor beta-3

Catalog No: #AP89042



Package Size: #AP89042-1 50ug #AP89042-2 100ug #AP89042-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	TGF-beta-3
Accession No.	Uniprot ID:P16047
Uniprot	P16047
GeneID	396438;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TGFB3

Note: This product is for in vitro research use only