

Tenascin-N

Catalog No: #AP89046



Package Size: #AP89046-1 50ug #AP89046-2 100ug #AP89046-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenascin-N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TN-N,Tenascin-W,TN-W,Tnw
Accession No.	Uniprot ID:Q80Z71
Uniprot	Q80Z71
GeneID	329278;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnn

Note: This product is for in vitro research use only