

## Telomerase reverse transcriptase

Catalog No: #AP89047



Package Size: #AP89047-1 50ug #AP89047-2 100ug #AP89047-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Telomerase reverse transcriptase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Telomerase catalytic subunit
Accession No.	Uniprot ID:Q6A548
Uniprot	Q6A548
GeneID	403412;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TERT

Note: This product is for in vitro research use only