

Vascular endothelial growth factor receptor 1

Catalog No: #AP89049



Package Size: #AP89049-1 50ug #AP89049-2 100ug #AP89049-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	VEGFR-1, VEGFR1
Accession No.	Uniprot ID:Q8QHL3
Uniprot	Q8QHL3
GeneID	374100;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FLT1

Note: This product is for in vitro research use only