

Angiopoietin-1 receptor

Catalog No: #AP89051



Package Size: #AP89051-1 50ug #AP89051-2 100ug #AP89051-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endothelial tyrosine kinase, HYK, STK1, Tunica interna endothelial cell kinase, Tyrosine kinase with Ig and EGF homology domains-2, Tyrosine-protein kinase receptor TEK, Tyrosine-protein kinase receptor TIE-2, mTIE2, p140 TEK, CD202b, Hyk, Tie-2, Tie2
Accession No.	Uniprot ID:Q02858
Uniprot	Q02858
GeneID	21687;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: Tek

Note: This product is for in vitro research use only