

Tissue factor pathway inhibitor

Catalog No: #AP89055



Package Size: #AP89055-1 50ug #AP89055-2 100ug #AP89055-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tissue factor pathway inhibitor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TFPI, Extrinsic pathway inhibitor, EPI, Lipoprotein-associated coagulation inhibitor, LACI
Accession No.	Uniprot ID:Q02445
Uniprot	Q02445
GeneID	29436;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tfp1

Note: This product is for in vitro research use only