

Thymic stromal lymphopoietin

Catalog No: #AP89059



Package Size: #AP89059-1 50ug #AP89059-2 100ug #AP89059-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thymic stromal lymphopoietin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thymic stroma-derived lymphopoietin
Accession No.	Uniprot ID:Q9JIE6
Uniprot	Q9JIE6
GeneID	53603;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tslp

Note: This product is for in vitro research use only