

Toll-like receptor 1

Catalog No: #AP89061



Package Size: #AP89061-1 50ug #AP89061-2 100ug #AP89061-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toll-like receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Toll/interleukin-1 receptor-like protein,TIL,CD281,KIAA0012
Accession No.	Uniprot ID:Q15399
Uniprot	Q15399
GeneID	7096;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TLR1

Note: This product is for in vitro research use only