

TNF receptor-associated factor 6

Catalog No: #AP89062



Package Size: #AP89062-1 50ug #AP89062-2 100ug #AP89062-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF receptor-associated factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	E3 ubiquitin-protein ligase TRAF6,RING-type E3 ubiquitin transferase TRAF6
Accession No.	Uniprot ID:Q3ZCC3
Uniprot	Q3ZCC3
GeneID	539124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TRAF6

Note: This product is for in vitro research use only