

TNF receptor-associated factor 6

Catalog No: #AP89063



Package Size: #AP89063-1 50ug #AP89063-2 100ug #AP89063-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF receptor-associated factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	E3 ubiquitin-protein ligase TRAF6, Interleukin-1 signal transducer, RING finger protein 85, RING-type E3 ubiquitin transferase TRAF6, RNF85
Accession No.	Uniprot ID:Q9Y4K3
Uniprot	Q9Y4K3
GeneID	7189;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TRAF6

Note: This product is for in vitro research use only