

## Ubiquitin D

Catalog No: #AP89065



Package Size: #AP89065-1 50ug #AP89065-2 100ug #AP89065-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubiquitin D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Diubiquitin,Ubiquitin-like protein FAT10,FAT10
Accession No.	Uniprot ID:O15205
Uniprot	O15205
GeneID	10537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:UBD

Note: This product is for in vitro research use only