

Ubiquitin D

Catalog No: #AP89066



Package Size: #AP89066-1 50ug #AP89066-2 100ug #AP89066-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Diubiquitin,Ubiquitin-like protein FAT10,Fat10
Accession No.	Uniprot ID:Q921A3
Uniprot	Q921A3
GeneID	29168;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ubd

Note: This product is for in vitro research use only