

Ubiquitin-protein ligase E3A

Catalog No: #AP89067



Package Size: #AP89067-1 50ug #AP89067-2 100ug #AP89067-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-protein ligase E3A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HECT-type ubiquitin transferase E3A,Oncogenic protein-associated protein E6-AP
Accession No.	Uniprot ID:O08759
Uniprot	O08759
GeneID	22215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ube3a

Note: This product is for in vitro research use only