

Tryptophan 5-hydroxylase 2

Catalog No: #AP89071



Package Size: #AP89071-1 50ug #AP89071-2 100ug #AP89071-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptophan 5-hydroxylase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neuronal tryptophan hydroxylase, Tryptophan 5-monoxygenase 2, NTPH
Accession No.	Uniprot ID:Q8IWU9
Uniprot	Q8IWU9
GeneID	121278;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TPH2

Note: This product is for in vitro research use only