

## L-dopachrome tautomerase

Catalog No: #AP89072



Package Size: #AP89072-1 50ug #AP89072-2 100ug #AP89072-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-dopachrome tautomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DCT,DOPAchrome conversion factor,DOPAchrome isomerase,DOPAchrome oxidoreductase,L-dopachrome Delta-isomerase,SLATY locus protein,Tyrosinase-related protein 2,TRP-2,Tyrp-2,Tyrp2
Accession No.	Uniprot ID:P29812
Uniprot	P29812
GeneID	13190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dct

Note: This product is for in vitro research use only