

Thyrotropin subunit beta

Catalog No: #AP89073



Package Size: #AP89073-1 50ug #AP89073-2 100ug #AP89073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Thyrotropin subunit beta |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Thyroid-stimulating hormone subunit beta, TSH-B, TSH-beta, Thyrotropin beta chain |
| Accession No. | Uniprot ID:P54828 |
| Uniprot | P54828 |
| GeneID | 403973; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TSHB

Note: This product is for in vitro research use only