

# X-box-binding protein 1

Catalog No: #AP89074



Package Size: #AP89074-1 50ug #AP89074-2 100ug #AP89074-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	X-box-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hepatocarcinogenesis-related transcription factor,HTF,Htf
Accession No.	Uniprot ID:Q9R1S4
Uniprot	Q9R1S4
GeneID	289754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Xbp1

Note: This product is for in vitro research use only