

Lymphotactin

Catalog No: #AP89075



Package Size: #AP89075-1 50ug #AP89075-2 100ug #AP89075-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphotactin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C motif chemokine 1,Cytokine SCM-1,Small-inducible cytokine C1,Ltn,Scyc1
Accession No.	Uniprot ID:P51672
Uniprot	P51672
GeneID	171371;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Xcl1

Note: This product is for in vitro research use only