

Synphilin-1

Catalog No: #AP89077



Package Size: #AP89077-1 50ug #AP89077-2 100ug #AP89077-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Synphilin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-synuclein-interacting protein
Accession No.	Uniprot ID:Q99ME3
Uniprot	Q99ME3
GeneID	67847;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sncaip

Note: This product is for in vitro research use only